
Your Cystic Fibrosis Application Design Document



Name: Lorcan Bermingham Student Number: C00196903

November 6, 2017 Tutor: Dr Joseph Keogh

Contents


System Architecture	4
Database Layout	5
Patient Detail Use Case's	6
Login	6
Medical Profile	7
Register	8
Food Diary	8
Activity	10
Spirometer	11
Logout	12
Doctor Detail Use Case's	13
Login	13
Manage Patients	13
Logout	14
System Sequence Diagram's	15
Patient Sequence Diagrams	15
Register	15
Login	16
Food Diary	17
Activity	18
Spirometer	19
Doctor Sequence Diagrams	20
Login	20
Manage Patient's	21
Screens	22
Splash Screen	22
Login	22
Register/Information	23
Home/Food Diary	23
Lunch	24
	24
Activity	25

Your Cystic Fibrosis Application | Name: Lorcan Bermingham | 17/04/2018

Moderate Activity	25
Spirometer	26

System Architecture

Database



Your CF Application

Database Layout

```
Userjson

1 {
2  "User" : {
3   "0868663769" : {
4     "DOB" : "24/01/89",
5     "EIR" : "RW109",
6     "Email" : "brennanSinead@gmail.com",
7     "Firstname" : "Sinead",
8     "Lastname" : "Brennan",
9     "Password" : "password"
10     },
11     "0868889991" : {
12     "DOB" : "28/09/88",
13     "EIR" : "RY809",
14     "Email" : "healylorcan@gmail.com",
15     "Firstname" : "Lorcan",
16     "Lastname" : "Bermingham",
17     "Password" : "password1"
18     }
19     }
20  }
```

Patient Detail Use Case's

Login

Name	Login
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when the patient opens the application Patient enters mobile number and password Patient selects login The application calls on firebase to check the user and returns the object Use case ends when patient is logged in
Alternatives	 3a password is incorrect Application display an error warning Return to step 2 3b phone number is incorrect Application displays error warning Return to step 2

Medical Profile

Name	Medical Profile
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when the patient opens the application Patient navigates to the medical profile tab A form loads for the patient to fill out Patient fill out the form The use case ends when the patient has finished filling out the form
Alternatives	 3a the form doesn't load Application display an error warning Return to step 3 3b the patient can't navigate to the medical profile tab Application displays error warning Return to step 2

Register

Name	Register
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when patient opens the application Patient selects register which then displays an information page Patient is then redirected to the register page Patient is prompted to enter account details Entered details are checked against users in firebase Once the details have been checked the patient is redirected to the login page The use case ends when the patient has logged in
Alternatives	 4a phone number already exists The application displays the number already exists Return to step 4 4b characters have been entered instead of numbers Application displays numbers only Return to step 4

Food Diary

Name	Food Diary
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when a patient has logged in Firebase loads stored image urls and names for the food diary Patient is brought to the home page which doubles as the food diary Patient can select breakfast, lunch, dinner or snacks Patient selects one of the abovementioned options Firebase loads the names and images from the database depending on which option was selected. Patient selects a food item from the selected food group. They are then met with a short description, there calorie counter for the day, the quantity of items they wish to choose and the specific calorie count of the food they have selected Once the item and the amount have been selected the application will automatically subtract that amount from there predefined amount set to the patient and update that amount in firebase along with the food details. This use case ends when the user logs out
Alternatives	 2a the images don't load when the user logs in The application displays empty cards only with the name The application tries to reload the images 4a the application fails to load the selected food area The patient selects a food group and the application fails to load the selected food group The application refreshes and tries again

Activity

Name	Activity
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when a patient has logged in Patient selects Activity from the menu Firebase loads stored image urls and names for the activity Patient is brought to the activity page Patient can select Light, Moderate, or Vigorous activities Patient selects one of the abovementioned options Firebase loads the names and images from the database depending on which option was selected. Patient selects an activity from the chosen activity group. The Patient is met with a short description and can select how long they want to do the chosen activity. Once the activity and the length have been selected the application will automatically multiple the length of time by the predefined number set to that activity and update that score in firebase along with the activity. The use case ends when the patient leaves the activity menu
Alternatives	 3a the images don't load when the user selects activity The application displays empty cards only with the name The application tries to reload the images 9a the application fails to load the selected activity area The patient selects a activity group and the application fails to load the selected activity group The application refreshes and tries again

9b the application makes the wrong calculation The application makes an error when calculating the activity's score The application flags the error and try's again.
and try's again

Spirometer

Name	spirometer
Actors	Patient, Firebase, Spirometer
Main Success Scenario	 This use case starts when the patient receives a push notification to do their spirometer The patient navigates to the spirometer tab in the menu The patient can optionally choose to follow a guide explaining how to use the spirometer The patient is redirected to spirometer page. The patient starts the timer and enters the result of each of the three tests they have completed. The patient stops the timer and clicks the submit button The application takes the three input values converts it from a String into a INT and the values are added together and divided by three to get the average of the three attempts. The new average and the time are then sent to the database This use case ends when the patient closes the spirometer tab
Alternatives	1a the patient doesn't receive a push notification The patient isn't informed to do the spirometer test The system notices that the test has been completed and tries

to resend the push notification • 5a the timer doesn't start The patient clicks the timer to
start and it falls to start
The system refreshes the page

Logout

Name	Logout
Actors	Patient, Firebase
Main Success Scenario	 The use case starts when the patient is logged in The patient selects the logout tab The use case ends when the patient is logged out
Alternatives	2a the logout tab doesn't work The patient selects the logout tab and it doesn't work Return to step 2

Doctor Detail Use Case's

Login

Name	Login
Actors	Patient, Firebase
Main Success Scenario	 Use case begins when the doctor opens the application Doctor enters username and password Doctor selects login The application calls on firebase to check the user and returns the object Use case ends when doctor is logged in
Alternatives	 3a password is incorrect Application display an error warning Return to step 2 3b phone number is incorrect Application displays error warning Return to step 2

Manage Patients

Name	Manage Patients
Actors	Doctor, Patient, Firebase
Main Success Scenario	 This use case begins when the doctor has logged in The doctor selects the manage patient tab The doctor selects the patient they wish to monitor The doctor then can review the selected patient's data

Alternatives	 1a the doctor cannot login The doctor receives an error when trying to login Return to step 1 3a the patients are not displayed The doctor selects the manage patient tab and noting is displayed Return to step 2


Logout

Name	Logout
Actors	Doctor, Firebase
Main Success Scenario	 The use case starts when the patient is logged in The patient selects the logout tab The use case ends when the patient is logged out
Alternatives	2a the logout tab doesn't work The patient selects the logout tab and it doesn't work Return to step 2


System Sequence Diagram's

Patient Sequence Diagrams


Register




Login




Food Diary



Activity




Spirometer




Doctor Sequence Diagrams

Login



Manage Patient's



Screens


Splash Screen



Login



Register/Information



Home/Food Diary




Lunch



Activity



Moderate Activity



Spirometer

